

Unique | Mine

U-MINE K



U-MINE

Beauty & Pain Care Combined in One

The World's first Microcurrent Amplification System



For Scalp, Face & Body Care
With Special Individual Accessories

“The World’s First All-in-One Beauty
Device”






U-MINE K
Material
& Functions

Dream Metal [Titanium]


Titanium is sturdy and has high strength, so it is used in aerospace and automotive industries. It is used as a raw material for structures and has excellent corrosion resistance, so is used in the chemistry industry.

In addition, it is widely used for shipbuilding, nuclear power, and industrial purposes to automobiles, golf clubs, mountaineering equipment, missiles, military vehicle armor, spacecraft, marine plants, seawater desalination, submarines, implants, artificial joints, artificial heart tubes, Kitchenware, and cosmetics.

U-Mine K is made from titanium, a raw material that is harmless to the human body!



Micro-current Flowing Through the Body



Cell membrane ion permeability increased by 40%

Helps tissue repair and regeneration of damaged cells

Collagen & elastin raw

Improves skin and restores elasticity, making skin younger

Increased ATP (adenosine triphosphate)

ATP, the source of energy for all cells, increases by 500% and it is excellent for removing waste products.

What if ATP is lacking?

Cell Necrosis and Aging

Skin aging, wrinkles, age spots, various inflammations

Cell Membrane Degeneration


Dementia, Alzheimer's, Parkinson's

Decreased Protein Synthesis

Immunity, Diabetes, Cardiovascular Disease

Decreased Biological Metabolism

Obesity, Hypothyroidism, Cancer, Autoimmune Disease and other
Chronic Diseases



Increase ATP by 500% with U-Mine K

ATP uses the energy contained in food to help cells complete their functions. It is the source of energy within cells that converts all necessary energy.

This molecule serves to transport energy.


It is present in cells of all body parts such as muscles, skin, brain, etc. ATP is responsible for promoting cellular energy production processes.



About 37
Trillion



40g per
minute ATP
consumption



2,400g per
hour ATP
consumption



About 58kg
ATP
consumption


ATP deficiency causes wrinkles and dryness on the skin when proteins stored in muscles are used to generate ATP.

 **U-Mine K can increase ATP by 500%!**

U-MINE K
Microcurrent
&
Solution Absorption



Solution Absorption




Microcurrent works on the principle that same poles push each other and different poles pull each other.

A method is used to ionize cosmetics into cathodes and anodes, and then naturally absorb them into the skin.

U-MINE K

Microcurrent & Pain Care



Pain Care

Diseases occur due to abnormalities in the potentials formed in the inner and outer membranes of cells in the body.



Paper on the Effects of Micro-Current on the Body



The Effects of Micro-Current on knee pain and arthritis



The Effects of Micro-Current on Stimulating Children's Growth Points



The Effects of Micro-Current on Stress Hormones



The Effects of Micro-Current on Chronic Pain

Articles on MicroCurrent Effects

chosun.com 사회

지면보기 구독신청 로그인 회원로그아웃

뉴스룸 오마이뉴스 정치 사회 경제 국제 스포츠 연예 한국 문화 라이프 스타일 해외 | 최정답뉴스

사회 달반 > 3개월 전기자극으로... 15년만에 깨어난 식물인간

조선일보 | 이명만 과학전문기자

14

0

711

입력 2017.09.26 01:04 | 수정 2017.09.26 09:21

[충남구민 원두 몸매 '전지약' 실효... 2개월만 움직임]

구글애드워드스 검색결과
 및 FDA 비안지오용 허가




연구진은 환자의 뇌를 안쪽에 있는 미주신경(迷走神経)에 전선을 얹고 전기자극을 줬다. 미주신경은 뇌와 인체의 모든 장기 사이를 오가며 신경 신호를 전달하는 통로이다. 통신망이 무너지면 사회에 불만이 발생하는 것처럼, 미주신경에 이상이 생기면 병이 난다. 전자약은 마치 통신망의 잡음을 제거하듯 인위적인 전기자극으로 잘못된 신경 신호를 교정해 치료 효과를 낸다.

전자약으로 식물인간 치료 방법

뇌 안쪽에 있는 미주신경에 전기 자극을 주면 뇌에서 운동·감각·의식을 담당하는 영역의 활동이 증가하면서 식물인간이 의식을 회복함

전기 자극 전후의 뇌 혈류량 변화



뇌에서 운동·감각·의식을 담당하는 두정엽과 시상 등에서 혈류량이 증가(붉은색) 현상을 얻을 수 있다.





U-MINE K Accessories

[Face & Scalp] Devices



✓ Face Lifting Effect
with Micro-Current & Sonic Vibration

✓ Skin Elasticity Effect by Promoting
Collagen

✓ Excellent Waste Discharge

✓ Helps Scalp Blood Circulation
to Prevent and Alleviate Hair Loss

✓ Prevention and Relief of Insomnia



[Body] Device



- ✓ Care for Pain in a Wide Area
- ✓ Damaged Nerve Recovery
- ✓ Visceral Fat Care & Diet Effects
- ✓ Sagging Muscle Lifting Effects
- ✓ Removal of Waste Product & Supply of Cell ATP

[Magic K] Device



- ✓ Powerful Pain and Inflammation Care
- ✓ Four Ergonomically Designed Feet
- ✓ Visceral Fat Care & Diet Effects
- ✓ Body Slimming & Skin Lifting
- ✓ Care for Women's & Men's issues



[K-Bell] Device



- ✓ Excellent Effects on Pain Care with Micro-Current & Sonic Vibration
- ✓ Discharge of Waste through Suction System
- ✓ Powerful Muscle Relaxation Pain Reduction Effect
- ✓ Home Care for Small and Large Areas



U-MINE K SYSTEM

U-MINE K

U-MINE K Main Body

Micro-Current & Sonic Vibration
Control System

Charge Indication
LED



Handle

Wired and Wireless Use Possible
Type C Charger

Current Power (1~6 level)


Power Button
Current & Vibration Control

Sonic Vibration Power



U-MINE KOREA's Special Technical Know-How

U-MINE K Accessories




Equipped with
Semiconductor
Current Conversion System
Current to MicroCurrent
Splitting System

Screw Connection
Easy to replace
accessories with
easy assembly

Accessory Function
Special Functions given to Each
Accessory

U-MINE K

Current Level 1·2 Effects



Main Effects


Face K Solution Absorption • Lifting Effect • Elasticity Effect

Scalp K Solution Absorption • Lifting Effect • Elasticity Effect •
Hair Loss Prevention • Scalp Circulation Improvement •
Dandruff Care

Body K & Pain Relief • Muscle Care
Magic K

U-MINE K

Current Level 3·4 Effects



Main Effects


Face K Solution Absorption • Lifting Effect • Elasticity Effect •
Wrinkle Improvement • Slimming Effect

Scalp K Hair Loss Prevention • Insomnia Improvement • Stress
Care

Body K &
Magic K Pain Relief • Muscle Care

U-MINE K

Current Level 5·6 Effects



Main Effects

Face K	Care for firm facial muscles Improvement of facial nerve paralysis Powerful slimming Powerful lifting Double chin improvement
Scalp K	Improvement of scalp problems such as nerve disease, inflammation Hair loss prevention Insomnia improvement Nerve recovery
Body K & Magic K	Rapid pain improvement Inflammation relief Promotion of healing process Improvement of areas with reduced sensation Accelerating the ligament healing progress

U-MINE K

Sonic Vibration Principle



1 & 2 Level

Decomposing waste

Reducing pores

Improving elasticity

Improving nerve stimulation

3 & 4 Level

Powerful slimming effect

Relaxing tight muscle

Bone correction

Powerful pain care

Powerful fat breakdown

Care Program

U-MINE K **Beauty Care Program**



Solving fundamental Scalp Problems

Scalp Lifting Care Program



Scalp-Accessory

Scalp lifting care for the fundamental lifting of the face prevents hair loss and induces comfortable sleep through anti-stress.



Scalp Lifting Care



Anti-Stress Care



Hair Loss Prevention
Insomnia Care

Face Care to Make your face look Smaller and Younger

Face Care Program



Face-Accessory

Cosmetics are absorbed deeply into the skin!
By caring the face line in detail,
you can fill the volume in the skin and
experience a lifting effect.



Slimming & Lifting Care
Cheekbone Reduction Care



Problem Skin Care
Solution Penetration Care



Wrinkle Care
Double Chin Care
Acne Pigmentation Care

For Soft Shoulder • Back • Neck

Back & Neck Care Program



Body-Accessory

Quickly Care for Hard and Sore
Shoulders and Neck!

Lightness and Freshness Experienced
after Care!



Muscle Relaxation
Shoulder and Neck Pain
Care



Blood Circulation &
Edema Care
Slim & Lifting Care



Inch Change Care
Muscle Relocation Care

Elasticity • Slim • Lifting • Lower Body Circulation Care

Lower Body Care Program




Elastic Lower Body Care by Lifting
Muscles

Slim • Lifting & Elasticity Care



Lower Body Blood
Circulation Care

Lower Body Swelling Care



Hip-Up Elasticity Care
Cellulite Tube



Slimming Care
Knee Care

Lipolysis Abdominal Lifting
Abdominal Care



**Rapid Fat Breakdown &
Reduction in Abdominal Size
Increase Breast Volume and Size**



Fat Breakdown Care
Abdominal Size Reduction



Abdominal Elasticity
Improvement Care
Chest Volume Up



Uterus Care
Surgery Scar Care

Certificates

방송통신기자재등의 적합등록 필증 Registration of Broadcasting and Communication Equipment	
상호 또는 명칭 Trade Name or Register	주미전자 유아정보회사
기자재명칭(제품명칭) Equipment Name	U-MINE K
기자재유형(기기종류) Equipment code	MS01 / LARS
기종도명칭 (Additional Equipment code)	
기종도번호 Main Model Number	UM-K20
과형도명칭 Series Model Number	
등록번호 Registration No.	K.R.11M-A-K-020
제조국(제조국가) Manufacture(Country of Origin)	주미전자 유아정보회사(중국)
등록일자 Date of Registration	2023-09-23
기타 (Others)	

본 기자재는 「방송법」 제18조제2항에 따라 등록의무를 증명합니다.
It is verified that foregoing equipment has been registered under the Clause 5, Article 18-2 of Radio Waves Act.

2023년 9월 23일(Monday) 22(일)까지

국립중앙연구원
Director General of National Radio Research Agency

KC Certificate

CE **DOC**
Declaration Of Conformity

We | U-MINERDIA CO., LTD.
121, Seomun-ro 389beon-gil, Gyeonggi-do,
Republic of Korea

DECLARE ON OUR SOLE RESPONSIBILITY THAT THE PRODUCT

Product name : U-MINE K
Model name : UM-K20

This CE Declaration of Conformity is following the provisions of RoHS (EU) 2011/65/EU and 2015/863/EU.
It is confirmed that a sample of the product has been tested and found to comply with below standards:

RoHS: RE 63000-2016, RE 63201-0-0-2013,
RE 63201-0-2013+AM01, 2017-030,
RE 63201-0-2013, RE 63201-0-2016,
RE 63201-0-1-2016, RE 63201-0-0-2017,
RE 63201-0-2017

(Detailed specification of the tested product shown in the following test report(s))
Model Test Report : N32001211206-0104

Date of issue : _____
Authorized By : _____

U-MINERDIA CO., LTD.
121, Seomun-ro 389beon-gil, Gyeonggi-do,
Republic of Korea

CE Certificate

Page 2 of 30

Test Report No.: NICE0011206-0104

RoHS TEST REPORT
European Directive 2011/65/EU, 2015/863/EU
Evaluation of RoHS Requirements for Electrical and Electronic Equipment

Trade mark:

Type of category: Category 04 Consumer equipment

Directive: 2011/65/EU, 2015/863/EU

Test standards: RE 63000-2016, RE 63201-0-1-2013, RE 63201-0-0-2013 + AM01, 2017-030, RE 63201-0-2013, RE 63201-0-2016, RE 63201-0-1-2016, RE 63201-0-0-2017, RE 63201-0-2017

Test result: The equipment which was evaluated has fulfilled with requirement of 2011/65/EU, 2015/863/EU Directive for the reasons:
RoHS: RE 63000-2016, RE 63201-0-1-2013, RE 63201-0-0-2013 + AM01, 2017-030, RE 63201-0-2013, RE 63201-0-2016, RE 63201-0-1-2016, RE 63201-0-0-2017, RE 63201-0-2017

Final test result : PASS

Approved by:

Natural Manager: Xiang Jun Chen

This publication may be reproduced in whole or in part for non-commercial purposes according to the of applicant. However, the user should refer to the report cover page without the absence of permission. Also, the publication requirement for the evaluation results of the second test item only any type of ESD, via full protocol, include assembly, adjustment or machine involving fourth test result. This publication results results only the tested item is compliant with RoHS equipment according to the evaluation procedures shown in this publication. The test report is not the unconditional result for formal laboratory accreditation systems, which signed the IAF-MLA. This report shall not be reproduced without the full, written approval of NICE Certification Co., Ltd.

Copyright © 2016 NICE Certification Co., Ltd. All rights reserved. NICE Certification Co., Ltd. is a member of the International Laboratory Accreditation Cooperation (ILAC) and the International Federation of Calibration and Testing Laboratories (IFCC).

ROHS Certificate

230010344759

检验报告
报告编号: CTL2308117041-Q

委托方: 厦门(天)生物科技有限公司
生产单位: 厦门(天)生物科技有限公司
样品名称: 锂电池保护
型号: UM-K20
检验类别: 委托检验

厦门(天)生物科技有限公司
2023.08.11

Chinese Certificate

U-MINE K Before & After





Before




After



Before



After



Before



After

PIC•COLLAGE



Before



After

PIC•COLLAGE





U-MINE K Models



U-MINE K Mini

1 Main Body, 1 K-Bell, 1 Scalp Accessory, 1 Charge
USD 1,500



U-MINE K Semi

1 Main Body, 1 K-Bell, 1 Face & 1 Scalp Accessory,
1 Body K, 1 Charge
USD 2,500



U-MINE K Special

2 Main Bodies, 2 K-Bell, 1 Magic K, 1 Body K

1 Face & 1 Scalp Accessory, 1 Charge

USD 4,000

Unique MINE

**Only One Brand
In the World**